

Fundamental Techniques and Measurements

- Mass Measurements
- Volume Measurements
- Preparation of a solution of known concentration
- UV-Visible Spectrophotometer

Electronic Balance

- What does an electronic balance measure?

- If you took an electronic balance with a capacity of 100 g to the moon what would its range be? _____

Why do good electronic balances have bubble level indicators?



Mass: Electronic Balance

- Accuracy
 - 4 to 6 significant digits
 - Calibration
 - Use known mass
 - Check weekly or when balance is moved
- Sources of error
 - Balance must be calibrated and maintained in same orientation in _____ field
 - hygroscopic chemicals: dry to constant mass first (will increase in mass rapidly as they reabsorb water on the balance!)
 - When preparing a solution of a given concentration it may be difficult to get the exact mass desired _____
 - evaporation of wet samples

Electronic Balance

| Model | Capacity | Resolution |
|---------|----------|------------|
| DI-100 | 100 g | 0.0001 g |
| DI-800 | 800 g | 0.01 g |
| DI-5000 | 5000 g | 0.1 g |



- For maximum accuracy use balance with _____ resolution possible!
- Don't forget to clean the balance if you spill any chemicals!!!!!!


Volume

- Volumetric flask
 - accuracy of _____/100 mL
 - Graduated cylinder
 - accuracy of _____/100 mL
 - Beaker
 - accuracy of _____/100 mL
 - Pipette
 - accuracy of \pm _____ for 100-1000 μ L
 - accuracy of \pm _____ for 10-100 μ L
- What will accuracy of solution be if you use pipette, volumetric flask, and electronic balance? _____
 What controls the accuracy? _____

Digital Pipettes

- Air displacement
- Do not directly contact fluid volume
 - avoids contamination of pipette
 - avoids sample carryover
- Require air tight connection between tip and body

Pipette Workings



Preparation of Solutions

- Example: Prepare 100 mL of a 30 mM solution of methylene blue.
- The molecular weight of methylene blue ($C_{16}H_{18}N_3SCl$) is 319.87 g.

$$CV = M$$

$$\underbrace{30 \times 10^{-3} \frac{\text{mole MB}}{\text{L}}}_{\text{concentration}} \cdot \underbrace{\frac{319.87 \text{ g MB}}{\text{mole MB}}}_{\text{conversion}} \cdot \underbrace{100 \times 10^{-3} \text{ L}}_{\text{volume}} = \underbrace{0.9596 \text{ g MB}}_{\text{mass}}$$

Preparation of Dilutions

- Prepare 100 mL of a 300 μM solution from the 30 mM solution
- Conservation of _____

$$M_{\text{dilute}} = M_{\text{concentrate}}$$

$$C_{\text{dilute}} V_{\text{dilute}} = C_{\text{concentrate}} V_{\text{concentrate}}$$

$$V_{\text{concentrate}} = \frac{C_{\text{dilute}} V_{\text{dilute}}}{C_{\text{concentrate}}}$$

$$V_{\text{concentrate}} = \frac{(300 \mu\text{M})(100 \text{ mL})}{(30 \text{ mM})}$$

Preparation of Solutions


- Fill volumetric flask half way with distilled water
- Add reagent (could be solid or liquid)
- Mix
- Fill volumetric flask to the line
- Mix
- Verify that volume didn't change (if necessary refill to line)



UV-Visible Spectrophotometer

- Theory
- Instrument
- Sample requirements
- Software

Light Attenuation by an Aqueous Solution



P is light intensity (photons/s)

$$\frac{dP}{dx} = -kP \quad k \propto C$$

$$\int_{P_0}^P \frac{dP}{P} = - \int_0^x k dx$$

$$\ln \left(\frac{P}{P_0} \right) = -kx$$

Theory: Light Attenuation = f(?)


- For a given excitation process, a molecule absorbs only one discrete amount of energy: expect very narrow absorption lines.
- Different vibrational and rotational states yield _____ absorption lines.
- First order decay with distance

$$\ln\left(\frac{P}{P_0}\right) = -kx \quad A = \log \frac{P_0}{P} = \epsilon bc$$

$$A = \epsilon bc$$

- P_0 - _____ light intensity
- P light intensity after passing through sample
- b - _____
- c - _____
- ϵ - _____ coefficient (function of wavelength and molecule)


Absorption Spectra



- Absorption Spectra for Methylene Blue


- Broad peaks
- Absorbs _____, looks _____

Instrument Light Path




Sample Requirements

- Sipper cell
 - peristaltic pump draws sample into sipper cell
 - requires a few mL to displace previous cell contents



Software



- Reference (single sample)
 - measures absorbance of sample cell and reference solution
 - usually distilled water or reagent blank
- Standards (multiple samples)
 - used to create a _____ curve
- Samples (multiple samples)
 - after sampling, standards can be used to estimate the concentration of samples

Absorbance Measurement Limitations


- P_o is a function of the ____.
- If absorbance is high what is P ? ____
- Suppose $A = 3$, what is P_o/P ? ____
- Suppose I create samples of higher and higher concentration. What will happen to the absorbance measurements?

$$A = \log \frac{P_o}{P} = \epsilon bc$$

There is a ____ (non zero) P that can be measured by an instrument.
 A ____ keep increasing! ____

Maximum Absorbance: P_o is measured as reference!

- Max absorbance $f()$
 - _____
 - _____
 - _____
 - _____
 - _____
- absorbance readings that exceed this value will not be used in analysis



Standards

- your name
- general description
- rinse time
- sample time
- sample concentrations
- select number of samples by moving this control

Samples

- enter sample descriptions here
- select number of samples by moving this control

