



## pH measurements

The probe, from chemistry to voltage to a number


## pH probe



## Glass membrane (Insulator)



## pH Measurements



- The porous frit provides electrical contact between the solution and the electrolyte
  - Must be in contact with the sample solution
  - Probe won't work well if frit is clogged (fouled)
- The voltage measurement requires a very high impedance circuit (high resistance) because a pH probe can't produce much current
- Gentle stirring keeps the solution next to the glass bulb from being depleted of protons

## Difficult Measurements?

- pH is difficult to measure in poorly buffered solutions
  - Distilled water
  - Rain
  - Between pKs of dilute buffers

## Nernst Equation: Voltage = f(pH)

$$E = E^0 + \frac{RT}{nF} \ln \left( \frac{[H^+]}{[H^{+0}]} \right) \quad [H^{+0}] \text{ Reference (known) } [H^+]$$

$$E = E^0 + \frac{RT \ln(10)}{nF} \log \left( \frac{[H^+]}{[H^{+0}]} \right) \quad E^0 \text{ Voltage at } [H^{+0}]$$

$$E = E^0 + \frac{RT \ln(10)}{nF} [pH^0 - pH] \quad F \text{ Faraday constant}$$

$$R \text{ Gas constant}$$

$$R = 8.314 \frac{J}{mol \cdot K}$$


$$F = 96500 \frac{Coulombs}{mol \cdot e^-}$$

$$n = 1 \text{ mol } e^-$$

$$pH = pH^0 - \frac{(E - E^0)nF}{RT \ln(10)}$$

## Nernst Equation

$$pH = pH^0 - \frac{(E - E^0)nF}{RT \ln(10)} \quad pH = pH^0 + \left( \frac{E^0}{T^0} - \frac{E}{T} \right) \frac{nF}{R \ln(10)}$$




$$R = 8.314 \frac{J}{mol \cdot K}$$

$$F = 96500 \frac{Coulombs}{mol \ e^-}$$

$$n = 1 \text{ mol } e^-$$

## Slope vs. Temperature




- Temperature compensation is important when temperature changes between samples!

## pH Calibration

- It would be possible to make a pH measurement without any calibration
  - Based on theoretical values
  - This is how the software recognizes buffers!
- Calibration accounts for non-ideal probe behavior (fouling) as well as electronic measurement errors
- It is important that buffers be used covering the range of pH measurements


## The Challenge

0 to 3 known buffers used as standards  
Temperature compensation optional



## Multiple Point Calibrations w/ Temperature Compensation

- How would you use this information to calculate pH?



Calculate E/T

Piecewise linear fit

$$pH = pH^0 + \left( \frac{E^0}{T^0} - \frac{E}{T} \right) \frac{nF}{R \ln(10)}$$